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# Processing instruction

## Important information

- Read through these instructions for use carefully before every use and keep them easily accessible for the user or the corresponding skilled personnel.

- Read through the warnings marked with this symbol carefully. Improper use of the products can result in serious injuries to the patient, user or third parties.

### Manufacturer

Elcon Medical Instruments GmbH  
Dr. Karl- Storz- Str.26  
D- 78532 Tuttlingen  
Germany

## Information

- The national regulations and standards for the reprocessing of medical products must be complied with.
- Mechanical cleaning is preferable since this leads to a more effective result.
- All products must be removed from the transport packaging before use; other packaging, protective caps and protective films are to be removed and disposed of in an environmentally-friendly manner.
- If it is not a one-piece instrument, the product must be disassembled to the furthest degree possible.

## Warnings and precautionary measures

- In the event of damage, even during use, the product must be immediately replaced. Any other use can lead to complications and/or endangerment of persons during the operation.
- For manual cleaning, never use metal brushes, metal sponges or abrasive cleaning agents.
- Defective products must have run through the entire reprocessing procedure before being sent back for repairs or claims.
- Never use cleaning solutions with bleach, such as sodium hypochlorite, since these can cause strong corrosion.
- The products must not be overstressed by bending, twisting or prying, since this can result in damage or breakage.

## Reprocessing restriction

- Frequent reprocessing has minor effects on the instruments. The product service life is limited by wear due to its intended use and damage to the instrument. The product may no longer be used under the following aspects (among others): Corrosion, damage, breaks, cracks, deformation, porosity, functional restrictions, products with illegible or missing label. Therefore, the corresponding information for checking the function must be complied with by the processor. Please properly dispose of or recycle the instruments after the product service life is over. The national regulations and disposal guidelines must be observed here.
- The service life of the products mainly depends how carefully the user handles them during use and reprocessing.

## Range of application

- The instruments may only be used for their intended purpose in the medical fields by correspondingly trained and qualified personnel. The treating physician / user is responsible for selecting the instruments for specific applications / surgical use, the appropriate training and information, and sufficient experience for working with the instruments.

## Safe handling and provision

### General

- All instruments from ELCON Medical Instruments GmbH are delivered unsterilized and must be cleaned, disinfected and sterilized by the user before use.
- The user should become familiar with the instrument and its functionality before use.
- The instructions for use provided by ELCON Medical Instruments GmbH must be read, complied with and kept in a safe place.
- Always use instruments as intended.
- Keep factory-new or unused instruments in a dry, clean, protected place.
- The instruments must be checked for their usability and functionality before every use.
- Never use damaged instruments. Do not make repairs yourself. Servicing and repairs should only be performed by appropriately qualified personnel.
- Only replace damaged components with original replacement parts.

### Reprocessing instructions

- Comply with the national legal regulations and standards for reprocessing.
- For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD cases or possible variants, comply with the respectively valid national ordinances regarding the reprocessing of the products.
- Purely manual cleaning – even using an ultrasonic bath – should only be applied if a mechanical method is not available, due to the much lower effectiveness and reproducibility.
- Note that the successful reprocessing of this medical product can only be ensured after the reprocessing process has been previously validated. The user/processor is responsible for this. Due to process tolerances, the manufacturer specifications serve only as a guide value for the reprocessing processes that exist at the processor site.
- For further current information about reprocessing, see also [www.a-k-i.org](http://www.a-k-i.org)

## Instrument reprocessing

### General information

Dried-on/fluxed operation residues can make cleaning more difficult or ineffective and can lead to corrosion in the case of stainless steel. Therefore, a time period between use and processing should not exceed 6 h, no fixing pre-cleaning temperatures > 45 °C should be applied, and no fixing disinfectants (active ingredient basis: aldehyde, alcohol) should be used.

Overdosed neutralizers or basic cleaners can lead to chemical attack and/or fading of the labelling on stainless steel.

In the case of stainless steel, chlorine or chloride-containing residues, e.g. contained within operation residues, tinctures, medicines, saline solutions, service water for cleaning, cleaning agents/disinfectants, can lead to corrosion damage (pitting corrosion, stress corrosion) and therefore to the destruction of the products. To remove, rinse adequately with deionized water and then dry.

Only process chemicals may be used which were recommended by the chemical manufacturer with regard to cleaning/disinfection effect as well as material compatibilities. All requirements for use, such as those concerning temperatures, concentrations, treatment times, etc. are to be strictly complied with. Otherwise, this can lead to the following problems:

Optical material changes, such as fading or colour changes on titanium or aluminium. In the case of aluminium, visible surface changes might already occur for a pH value >8 in the application solution, or material damage, such as corrosion, cracks, breaks, premature ageing or swelling. For further detailed information concerning hygienically safe reprocessing that protects and preserves the value of the material, see [www.a-k-i.org](http://www.a-k-i.org) under the red brochure publication "Instrument Reprocessing".

### Preparation before cleaning

- If it is not a one-piece instrument, the product must be disassembled to the furthest degree possible.
- Pre-clean instruments as needed.
- Reprocess instruments immediately after use.
- In the case of wet removal, use suitable cleaning agents/disinfectants. Before mechanical cleaning and disinfection, rinse thoroughly under running water.
- Rinse instruments with lumina 5x with a disposable syringe (minimum volume 50 ml).
- If necessary, clean product with ultrasound; see cleaning/disinfection.
- After the ultrasonic bath, rinse through the lumen for 10 seconds manually with a water pistol.

### Cleaning/disinfection

Damage to the product due to unsuitable cleaning agent/disinfectant and/or excessively high temperatures! Use cleaning agents and disinfectants which are approved for surgical steels according to the manufacturer instructions. Observe specifications concerning concentration, temperature and reaction time. Do not exceed the maximum permissible cleaning temperature of 55°C.

### Pre-cleaning

#### Perform ultrasonic cleaning:

- as effective mechanical support (recommended frequency 35 kHz) for manual cleaning/disinfection.
- for pre-cleaning products with dried-on residues before mechanical cleaning/disinfection, with 0.5 % neodisher Mediclean Forte for 10 min.
- as integrated mechanical support for mechanical cleaning/disinfection.
- for post-cleaning products with residues that were not removed after mechanical cleaning/disinfection.

Microsurgical products should be securely fixed to storage aids so they can be cleaned well in order to be able to continue with mechanical cleaning and disinfection

### Mechanical cleaning/disinfection with manual pre-cleaning

| Phase | Step  | T [°C/F]   | t [min] | Water quality | Chemical                          |
|-------|---|------------|---------|---------------|-----------------------------------|
| I     | Pre-rinsing   | cold       | 1       | -             | -                                 |
| II    | Cleaning  | 55°C ± 5°C | 5       | VE-W          | neodisher® Mediclean forte (0.5%) |
| (III) | (Neutralization can be omitted if neodisher® Mediclean forte is used) | -          | -       | -             | neodisher® Z                      |

| IV | Intermediate rinsing | RT         | 1  | VE-W | - |
|----|----------------------|------------|----|------|---|
| V  | Thermal disinfection | > 90°C     | 5  |      |   |
| VI | Drying               | 60°C ± 5°C | 30 | -    | - |

T-W: Drinking water

VE-W: Deionized (demineralized) water

- Place instruments on a sieve basket suitable for cleaning (avoid "rinsing shadows" - areas which could not be properly rinsed due to placement in cleaning container).
- It is recommended that deionized water be used in the cleaning step as well as in the final inspection. The neutralization step required for the classical alkaline cleaners can be omitted.
- The above-named application solution is to be completely rinsed off with water (preferably deionized).
- Inspect visible surfaces for residues after mechanical cleaning/disinfection.
- If necessary, repeat the cleaning process.

## Check, care and inspection

- Allow instruments to cool off to room temperature.
- Inspect the instruments after every cleaning and disinfection for: cleanliness, function and damage, e.g. insulation, loose, bent, broken, cracked, worn and broken-off parts.
- Immediately sort out damaged instruments.
- Maintenance means applying instrument oil or instrument milk (emulsion of white oil in water). Instruments with joints or closures (scissors, clamps, etc.) or with metallic sliding surfaces (rib shears, punches, etc.) must be treated with steam-sterilizable maintenance products with a paraffin oil base. The paraffin oil must comply with the respectively valid pharmacopoeia and be physiologically safe according to the German pharmacopoeia, 10. Issue (DAB 10) or the European pharmacopoeia (Ph. Eur.) or American pharmacopoeia, United States Pharmacopoeia (USP). The maintenance products prevent metal-on-metal friction and keep the instruments moving smoothly. Laser-labelled products can fade if treated with phosphoric acid- and hydrofluoric acid-containing basic cleaners. As a result, the coding function can be impaired or can be lost. Generally, surgical instruments must undergo constant maintenance before the function check. Maintenance products must work such that, even if they are constantly used, the joint parts do not stick together from the additive effect.

## Packaging

- If not done, the disassembled product must be reassembled.
- Protect instruments having a fine working end accordingly.
- Sort instruments into the associated storage location or to a suitable sieve basket.
- Package the sieve baskets appropriately for the sterilization method (make sure that the packaging prevents recontamination of the product between reprocessing and reuse).

## Sterilization

- Standardized packaging of instruments for sterilization in sterile bags in acc. with ISO 11607-1. A deviation from the suggested packaging must be validated.
- Make sure that the means of sterilization has access to all external and internal surfaces (e.g. by opening valves and taps).
- Validated sterilization method
- Steam sterilization in fractionated vacuum method (2 vacuum phases)
- Sterilization in fractionated vacuum method at temperatures of 132°C / holding time 5 min / drying time 10 min
- For simultaneous sterilization of multiple products in one steam sterilizer: Make sure that the maximum permissible loading of the steam sterilizer is not exceeded according to manufacturer specifications.

## Packaging

Use appropriate packaging for sterilization. (Either use transparent, sterile packaging or a suitable sterilization cassette (observe the instructions for use of the sterilization cassette manufacturer).

### Storage

Store sterilized instruments in a dry, clean, dust-free environment at moderate temperatures.

### Warranty

The instruments were made of high-quality material and underwent a quality check before delivery. If errors should nevertheless occur, please refer to the address given above. The warranty is void if repairs are carried out by companies not authorized for this by Elcon Medical Instruments GmbH.

### Servicing and repairs

Do not make any unauthorized repairs or modifications to the product. Only authorized personnel of the manufacturer is responsible and intended for this. If you should have any complaints, claims or information regarding our products, we ask that you contact us.

### Return transport

Defective or non-compliant products must have run through the entire reprocessing procedure before being sent back for repairs/service.

## Inspection instruction

Before every instrument use, it is to be examined for breaks, cracks, deformations and proper functioning. Worn, corroded, deformed, porous or otherwise damaged instruments must be sorted out.

- The following components must be checked especially carefully:
  - Notches
  - Blades
  - Tips
  - Joints
  - Handle springs

The stainless steels (non-tasting) used for instrument manufacture form specific passive layers as protective layers based on their alloy. These steels are only conditionally resistant to the attack of chloride ions and aggressive media and liquids. In addition to the efforts taken by the manufacturer in selecting the right materials and their diligent processing, the instruments must be professionally and continuously maintained and processed properly at the user site.

## Material resistance

When selecting the cleaning agent and disinfectant, make sure that the following ingredients are not contained in them:

- Organic, mineral and oxidizing acids (minimum permissible pH value is 5.5).
- Bases/strong bases (neutral/enzymatic (max. permissible pH value 8.5, mandatory for products made of aluminium or other alkal-sensitive materials; see the chapter "Special instructions") or alkaline cleaners (max. permissible pH value 11, mandatory for products with intended use in prion-critical areas, e.g. according to Attachment 7 of the KRINKO RKI BIArM recommendation for processing) recommended)
- Organic solvents (e.g., alcohols, ethers, ketones, benzenes)
- Oxidation agents (e.g., hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- Aromatichalogenic hydrocarbons

## Used symbols

|  |   |     |                                  |
|--|---|-----|----------------------------------|
|  | ATTENTION<br>(Attention, observe accompanying documents") | REF | Order number or catalogue number |
|  | Observe the instructions for use                          | LOT | Batch designation                |
|  | Manufacturer  |     | Attention: Non-sterile product   |
|  | CE identification   |     |                                  |

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