



Australian Register of Therapeutic Goods Certificate

Issued to

InstruMED

for approval to supply

InstruMED - Osteotome - Osteotome

ARTG Identifier	305737
ARTG Start Date	29/06/2018
Product Category	Medical Device Included Class 1
GMDN	12844
GMDN Term	Osteotome
Intended Purpose	A surgical instrument having a cutting edge that is bevelled on either side with sloping curves. It is used during an osteotomy (the surgical procedure to divide bone); this bone will typically be used as a transplant. The osteotome may also be used for shaping bone.

Manufacturer Details	Address	Certificate number(s)
Elcon Medical Instruments GmbH	Dr Karl-Storz-strasse 26A Tuttlingen , 78532 , Germany	

ARTG Standard Conditions

The above Medical Device Included Class 1 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Osteotome - Osteotome

Product Specific Conditions

No specific conditions have been recorded against this entry.