



Australian Register of Therapeutic Goods Certificate

Issued to

InstruMED

for approval to supply

InstruMED - Punch, surgical, bone - Punch, surgical, bone

ARTG Identifier	305848
ARTG Start Date	29/06/2018
Product Category	Medical Device Included Class 1
GMDN	35287
GMDN Term	Punch, surgical, bone
Intended Purpose	A surgical instrument used to punch holes in bone, typically for the purposes of biopsy procedures or for fixation procedures. The instrument can be ring handled, pistol grip or shaft-like in design. Ring handled or pistol grip designs feature extended shafts which terminate in a punch-like or mechanical jaw mechanism at the distal tip. The shaft-like instrument is a hollow tube, sharpened on one end, and is usually used with a trocar.

Manufacturer Details	Address	Certificate number(s)
Elcon Medical Instruments GmbH	Dr Karl-Storz-strasse 26A Tuttlingen , 78532 , Germany	

ARTG Standard Conditions

The above Medical Device Included Class 1 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Punch, surgical, bone - Punch, surgical, bone

Product Specific Conditions

No specific conditions have been recorded against this entry.