



Australian Register of Therapeutic Goods Certificate

Issued to

InstruMED

for approval to supply

InstruMED - Rasp, surgical, bone - Rasp, surgical, bone

ARTG Identifier	305852
ARTG Start Date	29/06/2018
Product Category	Medical Device Included Class 1
GMDN	35559
GMDN Term	Rasp, surgical, bone
Intended Purpose	A metal, file-like surgical instrument with a ridged surface of coarse grooves or teeth designed to be used in surgical procedures, e.g. orthopaedic surgery. It is used to scrape and smooth the surface of bone. A rasp differs from a file in that the grooves or teeth are much coarser so that they will not clog so easily (as would a file) and render the instrument unusable until cleaned of the debris.

Manufacturer Details	Address	Certificate number(s)
Elcon Medical Instruments GmbH	Dr Karl-Storz-strasse 26A Tuttlingen , 78532 , Germany	

ARTG Standard Conditions

The above Medical Device Included Class 1 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Rasp, surgical, bone - Rasp, surgical, bone

Product Specific Conditions

No specific conditions have been recorded against this entry.