



Australian Register of Therapeutic Goods Certificate

Issued to

InstruMED

for approval to supply

InstruMED - Wire holding/twisting pliers - Wire holding/twisting pliers

ARTG Identifier	305881
ARTG Start Date	29/06/2018
Product Category	Medical Device Included Class 1
GMDN	32874
GMDN Term	Wire holding/twisting pliers
Intended Purpose	A heavy-duty, surgical instrument with specially designed sturdy handles and blade jaws used to grip, tighten and/or twist wires that are being applied to the patient during a surgical intervention. It has a scissors-like design with curved handles and is made of high-grade stainless steel. It is available in various sizes and the jaws (typically short and broad with serrations) at the working end may be pivoted or have a ratio-lever exchange to provide greater gripping force. This is a reusable device.

Manufacturer Details	Address	Certificate number(s)
Elcon Medical Instruments GmbH	Dr Karl-Storz-strasse 26A Tuttlingen , 78532 , Germany	

ARTG Standard Conditions

The above Medical Device Included Class 1 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Wire holding/twisting pliers - Wire holding/twisting pliers

Product Specific Conditions

No specific conditions have been recorded against this entry.