

## Quality that you can track



Each instrument is laser marked with:

- > Manufacturer (ELCON)
- > Manufacturing origin (GERMANY)
- > Item Number
- > Material (Stainless)
- > Lot Number
- > UDI Data Matrix Code



The data matrix code is compliant and contains all information required in Europe (EN ISO 13485) and the USA (FDA).

The information is as follows:

UDI - DI: Unique Device Identification (DI = Device Identifier)

This is the GTIN (Global Trade Item Number) and it is for every instrument unique.

UDI - PI: Unique Device Identification (PI = Production Identifier)

This is the ELCON Lot Number

An industry-leading TRUMPF laser marking machine is used by the factory in Germany



### CONTACT



**08 8354 0525**

Call for more information



[sales@instrumed.com.au](mailto:sales@instrumed.com.au)



[www.instrumed.com.au](http://www.instrumed.com.au)